

Anti- *C. perfringens* collagenase antibody, mouse monoclonal (cp-02)

64-050 100 µg

Storage: -20°C.

Immunogen: Culture supernatant of *Clostridium perfringens*

Form: 1 mg/ml in PBS⁻ with 50% glycerol, filter sterilized.

Purity: Purified using Ab-Capture for IgM (ProteNova, Japan)

Isotype: mouse IgM

Reactivity: Reacts with collagenases of *Clostridium perfringens* and *C. histolyticum*

Applications:

1. Western blotting (1/500~1/1,000)
2. ELISA (assay dependent)

This antibody is useful for detecting food-poisoning *Clostridium* strains.

Background: *Clostridium perfringens* is one of the major causative agents of food poisoning. *C. perfringens* produces various gelatinolytic enzymes with molecular masses ranging from approximately 120 to approximately 60 kDa. A gelatinolytic enzyme is present in the largest quantity in the culture supernatant, and this enzyme is purified as collagenase. The collagenase of *Clostridium histolyticum* (68 kDa) is the best studied and characterized.

Data Link: UniProtKB: [P43153](#) ((COLA_CLOPE), [Q46085](#) ((COLH_HATHI)




Fig.1. Detection of collagenase of *C. perfringens* by Western blotting with monoclonal antibody (MAb cp-02).

1. Purified collagenase of *C. histolyticum*
2. Culture supernatant of *C. perfringens*.

The 68 kDa band in lane 2 is collagenase of *C. perfringens*.

The primary antibody was used at 1/1,000 dilution.





Fig.2. Titration of antibody reactivity of MAb (cp-02) by indirect ELISA using culture medium of *C. perfringens*.

The wells of plate were coated with culture medium of *C. perfringens* (100 μ l, 1 μ g/ml). After blocking with 5% skim milk, 100 μ l of antibody at the indicated dilutions was added to the each well. HRP-conjugated goat anti-mouse IgG, IgM and IgA (100 μ l, x2000 dilution) was added. Color was developed with orthophenylenediamine as substrate. Optical densities (OD) measured at 490nm.

Fig.3. Titration of collagenase in culture medium of *C. perfringens* by indirect ELISA using MAb (cp-02).

ELISA plate is coated with indicated amounts of the culture medium of *C. perfringens* per well. MAb (cp-02) was used at 1/500 dilution. ELISA was performed as in Fig.2.



Tale 1. Immunological reactivity of MAb (cp-02) with various food poisoning bacteria

	ELISA	WB
<i>Clostridium perfringens</i> (ATCC13124)	+	68K
<i>Bacillus cereus</i>	—	
<i>Staphylococcus aureus</i>	—	
<i>Campylobacter jejuni</i>	—	
<i>Salmonella Enteritidis</i>	—	
<i>Vibrio parahaemolyticus</i>	—	
<i>Escherichia coli</i> (ETEC)	—	
<i>E. coli</i> 0157:H7 (EHEC)	—	
Purified Collagenase (from <i>C.histoliticum</i>)	+	68K

Reference: There has been no publication using this antibody.

Please note: All products are FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES. NOT FOR MILITARY USE.